



Developers Interface Guide

For Motorized Zoom Systems and
BrightLight[®] Controllers

*Visual C++ and Visual Basic.net
Serial and USB Connection*



Developers Interface Guide

Table of Contents

1. Introduction	5
1.1. Other Documents to reference	5
1.2. Products Covered	5
2. Physical Connection Environments.....	5
2.1. RS232	5
2.2. USB 2.....	5
3. Runtime considerations	5
3.1. Dynamic Link Library (DLL)	5
4. Development Environments	6
4.1. Visual C++	6
4.2. Visual Basic.NET.....	6
5. Sample code naming convention.....	6
5.1. Serial Port Connection Layout	7
6. Application Basics	7
6.1. DLL linkage.....	7
6.1.1. Visual C++	7
6.1.2. Visual Basic.NET	8
6.2. Connection initialization.....	8
6.2.1. Serial interface Visual C++.....	8
6.2.2. USB interface Visual C++	9
6.2.3. Serial interface Visual Basic.NET	9
6.2.4. USB interface Visual Basic.NET.....	10
7. Function call definitions.....	11
7.1. Serial interface.....	11
7.1.1. SerConnectionConnect.....	11
7.1.2. SerConnectionEstablished.....	11
7.1.3. SerConnectionDisconnect	12

Developers Interface Guide

Motorized Zoom Systems and BrightLight[®] Controllers

Page 3

7.1.4.	SerConnectionDisconnectAll	12
7.1.5.	SerConnectionWrite	12
7.1.6.	SerConnectionRead	12
7.2.	USB interface	13
7.2.1.	USBFindUSBInterfaces.....	13
7.2.2.	USBConnectionConnect	13
7.2.3.	USBConnectionConnectSpec.....	13
7.2.4.	USBConnectionEstablished	14
7.2.5.	USBConnectionDisconnect.....	14
7.2.6.	USBConnectionDisconnectAll.....	14
7.2.7.	USBConnectionWrite	15
7.2.8.	USBConnectionRead.....	15
8.	Motor Controller Register Definitions	15
8.1.	Motor Controller Identification Registers	15
8.1.1.	REG_SYS_PRODUCTID	15
8.1.2.	REG_SYS_VERSIONHW	16
8.1.3.	REG_SYS_VERSIONDATE	16
8.1.4.	REG_SYS_VERSIONSW	16
8.2.	Motor Controller System Setup Registers	16
8.2.1.	REG_SETUP_ACCEL.....	17
8.2.2.	REG_SETUP_INITVELOCITY.....	17
8.2.3.	REG_SETUP_MAXVELOCITY	18
8.2.4.	REG_SETUP_REVBACKLASH	18
8.2.5.	REG_SETUP_FWDBACKLASH	18
8.2.6.	REG_SETUP_CONFIG	19
8.2.7.	REG_SETUP_WRITE	19
8.3.	Motor Controller Operation Registers	20
8.3.1.	REG_USER_TARGET	20
8.3.2.	REG_USER_INCREMENT	20
8.3.3.	REG_USER_CURRENT	20

Developers Interface Guide

Motorized Zoom Systems and BrightLight[®] Controllers	Page 4
8.3.4. REG_USER_LIMIT.....	21
8.3.5. REG_USER_STATUS	21
8.3.6. REG_SETUP_LIMIT.....	22
9. BrightLight Register Definitions	22
9.1. BrightLight Identification Registers	22
9.1.1. REG_SYS_PRODUCTID	22
9.1.2. REG_SYS_VERSIONHW	23
9.1.3. REG_SYS_VERSIONDATE	23
9.1.4. REG_SYS_VERSIONSW.....	23
9.2. BrightLight Operation Registers	23
9.2.1. REG_PWM_ABSOLUTE	23
9.2.2. REG_PWM_INCREMENT	24



Developers Interface Guide

1. Introduction

1.1. Other Documents to Reference

Motor Controller User's Manual- Windows.doc

1.2. Products Covered

- Navitar 2 Phase Stepper Controller
- Navitar 5 Phase Stepper Controller
- Navitar DC Encoded Controller
- Navitar Digital BrightLight® Controllers

2. Physical Connection Environments

2.1. RS232

The covered products can be connected to a PC via an RS232 interface. The supplied **navserAPI.dll** configures the serial port with the necessary parameters and provides the connection.

2.2. USB 2

The covered products can be connected to a PC via a USB 2 interface. The supplied **navusbAPI.dll** provides the connection.

3. Runtime Considerations

3.1. Dynamic Link Library (DLL)

The low level communication interface for both the RS232 and USB environments is handled with a **Dynamic Link Library (DLL)**. The proper DLL must reside in the runtime path. The recommended location is to place the DLL in the same directory as the compiled executable.

Developers Interface Guide

4. Development Environments

4.1. Visual C++

Microsoft Visual C++ version 6.0 with service pack 5 or higher.

4.2. Visual Basic.NET

Microsoft Visual Basic.NET 2003 or higher.

5. Sample Code Naming Convention

There are 16 sample programs. The Visual Basic programs are located under the directory samples VB. The Visual C++ programs are located under the directory VC. Under each directory are 4 serial interface samples and 4 USB samples. Within each group of 4 samples there are the following samples:

- 1) BrightLight controller only
- 2) Motor controller only
- 3) 1 Motor controller and 1 BrightLight controller
- 4) 2 Motors controller and 2 BrightLight controller

The directory and project names are specified as follows:

Environment	Interface	Number of Motor controllers	Number of BrightLight Controllers	Folder name
VB	USB	0	1	VB/ubm0b1
VB	USB	1	0	VB/ubm1b0
VB	USB	1	1	VB/ubm1b1
VB	USB	2	2	VB/ubm2b2
VB	serial	0	1	VB/sbm0b1
VB	serial	1	0	VB/sbm1b0
VB	serial	1	1	VB/sbm1b1
VB	serial	2	2	VB/sbm2b1
VC	USB	0	1	VC/ubm0b1
VC	USB	1	0	VC/ubm1b0
VC	USB	1	1	VC/ubm1b1
VC	USB	2	2	VC/ubm2b2
VC	serial	0	1	VC/sbm0b1
VC	serial	1	0	VC/sbm1b0
VC	serial	1	1	VC/sbm1b1
VC	serial	2	2	VC/sbm2b2

Developers Interface Guide

5.1. Serial Port Connection Layout

The serial port samples use the following COM numbers:


Number of Motor Controller s	Number of BrightLight Controller s	COM1	COM2	COM3	COM4
0	1	BrightLight			
1	0	Motor			
1	1	Motor	BrightLight		
2	2	Motor	Motor	BrightLight	BrightLight

6. Application Basics

6.1. DLL linkage

6.1.1. Visual C++


The application should include the header file navserAPI.h or navusbAPI.h. Under Project/Settings Link tab the library file (**navserAPI.lib** or **navusbAPI.lib**) must be included.



Developers Interface Guide

6.1.2. Visual Basic.NET

The DLL interface description (**navserapi.vb** or **navusbapi.vb**) must be included in the project file.



6.2. Connection Initialization

6.2.1. Serial interface Visual C++

In order to establish a connection to a device the following code is used:

```
serHandle1 = SerConnectionConnect (COM1, DEF_MOTOR_CONTROLLER) ;  
serHandle2 = SerConnectionConnect (COM2, DEF_BRIGHTLIGHT) ;
```

The returned serHandle will be used in future calls to the device.

See the definition of SerConnectionConnect for more details.

Also examine the sample code near the end of the **OnInitDialog** function.

6.2.2.USB interface Visual C++

In order to establish a connection to a device the following line of code is used:

```
USBFindUSBInterfaces();  
  
usbHandle1 = USBConnectionConnect(1, DEF_MOTOR_CONTROLLER);  
usbHandle2 = USBConnectionConnect(1, DEF_BRIGHTLIGHT);
```

The returned usbHandle will be used in future calls to the device.

See the definition of USBConnnectionConnect for more details.

The call to USBFindUSBInterfaces must be called first. Then all the desired devices can be opened. The 'handles should be check for validity. If any handles are invalid, such as requested device not being connected or powered up the situation should be corrected and the sequence started again with the USBFindUSBInterfaces call.

If the application is running two different kinds of motor controllers the following is a more defined method of connecting to the correct motors:

```
USBFindUSBInterfaces();  
  
usbHandle1 =  
USBConnectionConnectSpec(1, DEF_MOTOR_CONTROLLER, DEF_MICROMO2PHASE);  
  
usbHandle2 =  
USBConnectionConnectSpec(1, DEF_MOTOR_CONTROLLER, DEF_VEXTA5PHASE);  
  
usbHandle3 = USBConnectionConnect(1, DEF_BRIGHTLIGHT);
```

Examine the sample code near the end of the **OnInitDialog** function.

6.2.3.Serial interface Visual Basic.NET

In order to establish a connection to a device the following line of code is used:

```
serHandle1 = SerConnectionConnect(COM1, DEF_MOTOR_CONTROLLER)  
serHandle2 = SerConnectionConnect(COM2, DEF_BRIGHTLIGHT)
```

The returned serHandle will be used in future calls to the device.

See the definition of SerConnnectionConnect for more details.

Also examine the sample code near the beginning of the:

```
#Region " Windows Form Designer generated code "  
after the InitializeComponent() call.
```

Developers Interface Guide

6.2.4.USB interface Visual Basic.NET

In order to establish a connection to a device the following line of code is used:

```
USBFindUSBInterfaces()  
  
usbHandle1 = USBConnectionConnect(1,DEF_MOTOR_CONTROLLER)  
usbHandle2 = USBConnectionConnect(1,DEF_BRIGHTLIGHT)
```

The returned usbHandle will be used in future calls to the device.

See the definition of USBConnnectionConnect for more details.

The call to USBFindUSBInterfaces must be called first. Then all the desired devices can be opened. The handles should be check for validity. If any handles are invalid, such as requested device not being connected or powered up the situation should be corrected and the sequence started again with the USBFindUSBInterfaces call.

```
USBFindUSBInterfaces()  
  
' connect to the 1st MicroMo motor found  
usbHandle1 =  
USBConnectionConnectSpec(1, DEF_MOTOR_CONTROLLER, DEF_MICROMO2PHASE)  
  
' connect to the 1st VEXTA motor controller found  
usbHandle2 =  
USBConnectionConnectSpec(1, DEF_MOTOR_CONTROLLER, DEF_VEXTA5PHASE)  
  
' connect to the 1st BrightLight found  
usbHandle3 = USBConnectionConnect(1, DEF_BRIGHTLIGHT)
```

Also examine the sample code near the beginning of the:

```
#Region " Windows Form Designer generated code "  
after the InitializeComponent() call.
```

Developers Interface Guide

7. Function Call Definitions

7.1. Serial Interface

7.1.1. SerConnectionConnect

int SerConnectionConnect(int **port**, unsigned short **ProductID**);

port	Specifies which COM port to connect to.
ProductID	Specifies which type of controller to connect to.
	DEF_BRIGHTLIGHT
	DEF_MOTOR_CONTROLLER
returns	An id to be used as a handle to future calls.

This routine attempts to open a connection to a specified device on the specified port.

If 0 is returned it indicates the call failed. Causes could be that the port could not be opened or that the specified controller was not found on the port.

7.1.2. SerConnectionEstablished

int SerConnectionEstablished(int **id**, unsigned short **ProductID**);

id	Handle returned by SerConnectionConnect.
ProductID	Specifies which type of controller to check.
	DEF_BRIGHTLIGHT
	DEF_MOTOR_CONTROLLER
returns	Connection status.

This routine asks if a controller of a specified type is connected to the specified port. This returns TRUE if the specified controller is connected, otherwise it returns FALSE.

7.1.3. SerConnectionDisconnect

```
int SerConnectionDisconnect(int id);
```

id	Handle returned by SerConnectionConnect.
-----------	--

This routine disconnects the device on the specified handle.

7.1.4. SerConnectionDisconnectAll

```
int SerConnectionDisconnectAll(void);
```

This routine disconnects all devices opened with SerConnectionConnect.

7.1.5. SerConnectionWrite

```
SerConnectionWrite(int id, unsigned char reg, long *v);
```

id	Handle returned by SerConnectionConnect.
reg	Register number.
*v	Pointer to the value to be written to the register.

This call writes the value pointed to by **v** to the register **reg**.

See the section on Register Definitions for more details.

7.1.6. SerConnectionRead

```
int SerConnectionRead(int id, unsigned char reg, long *v);
```

id	Handle returned by SerConnectionConnect.
reg	Register number.
*v	Pointer to the value to be read from the register.

This routine reads a value to the location pointed to by **v** from register **reg**.

See the section on Register Definitions for more details.

Developers Interface Guide

7.2. USB Interface

7.2.1.USBFindUSBInterfaces

```
int USBFindUSBInterfaces(void);
```

This routine parses the USB chain looking for Navitar devices.

7.2.2.USBConnectionConnect

```
int USBConnectionConnect(int device,unsigned short ProductID);
```

device	Specifies which of the specified product on the USB chain to connect to.
ProductID	Specifies which type of controller to connect to.
	DEF_BRIGHTLIGHT
	DEF_MOTOR_CONTROLLER
returns	An id to be used as a handle to future calls.

This routine attempts to open a connection to a specified device on the USB chain.

If 0x100 is returned it indicates the call failed. Causes could be that specified controller was not found on the port.

Example:

To connect to the first motor controller on the USB chain:

```
handle = USBConnectionConnect(1,DEF_MOTOR_CONTROLLER);
```

7.2.3.USBConnectionConnectSpec

```
int USBConnectionConnectSpec(int device,unsigned short ProductID,unsigned short ProductType);
```

device	Specifies which of the specified product on the USB chain to connect to.
ProductID	Specifies which type of controller to connect to.
	DEF_BRIGHTLIGHT
	DEF_MOTOR_CONTROLLER
returns	An id to be used as a handle to future calls.

Developers Interface Guide

This routine attempts to open a connection to a specified device on the specified port.

If 0x100 is returned it indicates the call failed. Causes could be that specified controller was not found on the port.

Example:

To connect to the first motor controller on the USB chain:

```
handle = USBConnectionConnect(1, DEF_MOTOR_CONTROLLER, DEF_VEXTA5PHASE);
```

7.2.4.USBConnectionEstablished

int USBConnectionEstablished(int **device**, unsigned short **ProductID**);

device	Handle returned by USBConnectionConnect.
ProductID	Specifies which type of controller to check.
	DEF_BRIGHTLIGHT
	DEF_MOTOR_CONTROLLER
returns	connection status

This routine asks if a controller of a specified type is connected to the specified port.

This returns TRUE if the specified controller is connected, otherwise it returns FALSE.

7.2.5.USBConnectionDisconnect

int USBConnectionDisconnect(int **id**);

id	Handle returned by USBConnectionConnect
-----------	---

This routine disconnects the device on the specified handle.

7.2.6.USBConnectionDisconnectAll

int USBConnectionDisconnectAll(void);

This routine disconnects all devices opened with USBConnectionConnect.

Developers Interface Guide

7.2.7.USBConnectionWrite

USBConnectionWrite(int **device**, unsigned char **reg**, long ***v**);

device	Handle returned by USBConnectionConnect
reg	Register number
*v	Pointer to the value to be written to the register

This call writes the value pointed to by **v** to the register **reg**.

See the section on Register Definitions for more details.

7.2.8.USBConnectionRead

int USBConnectionRead(int **device**, unsigned char **reg**, long ***v**);

device	Handle returned by USBConnectionConnect
reg	Register number
*v	Pointer to the value to be read from the register

This routine reads a value to the location pointed to by **v** from register **reg**.

See the section on Register Definitions for more details.

8. Motor Controller Register Definitions

The section entitled “Motor Controller Operation Registers” is probably the most interesting.

8.1. Motor Controller Identification Registers

8.1.1.REG_SYS_PRODUCTID

REG_SYS_PRODUCTID(0x01) read only

return value:

either

DEF_BRIGHTLIGHT 0x4000

DEF_MOTOR_CONTROLLER 0x4001

8.1.2.REG_SYS_VERSIONHW

REG_SYS_VERSIONHW(0x02) read only

return value:

32 bit hardware version identifier.

8.1.3.REG_SYS_VERSIONDATE

REG_SYS_VERSIONDATE(0x03) read only

return value:

32 bit build date identifier.

0x00,YY,MM,DD

example:

0x00050402 would be 2005, April, 2

8.1.4.REG_SYS_VERSIONSW

REG_SYS_VERSIONSW(0x04) read only

return value:

32 bit software version identifier.

8.2. Motor Controller System Setup Registers

These system registers are intended to be system setup registers. Changing these may result in a system that does not perform properly. The default values for these registers are listed in the document: Motor Controller User's Manual-Windows.

8.2.1.REG_SETUP_ACCEL

REG_SETUP_ACCEL_1(0x15) sets motor 1 acceleration

REG_SETUP_ACCEL_2(0x25) sets motor 2 acceleration

For Two and Five phase motors, these 2 registers set the acceleration value. If this value is set too high it is possible that the motor may stall.

For DC Encoded motors, these 2 registers set the velocity gain component. The predominant effect of this is controlling the deceleration of the motor.

WARNING: Changing these values from the factory defaults may result in a system that does not perform correctly.

See the document **Motor Controller User's Manual-Windows.doc** for the proper factory defaults.

8.2.2.REG_SETUP_INITVELOCITY

REG_SETUP_INITVELOCITY_1(0x16) sets initial motor 1 velocity

REG_SETUP_INITVELOCITY_2(0x26) sets initial motor 2 velocity

For Two and Five phase motors, these 2 registers set the initial step rate. If this value is set too high it is possible that the motor may stall on startup.

For DC Encoded motors, these 2 registers set the error gain parameter.

WARNING: Changing these values from the factory defaults may result in a system that does not perform correctly.

See the document **Motor Controller User's Manual-Windows** for the proper factory defaults.

8.2.3.REG_SETUP_MAXVELOCITY

REG_SETUP_MAXVELOCITY_1(0x17) sets maximum motor 1 velocity

REG_SETUP_MAXVELOCITY_2(0x27) sets maximum motor 2 velocity

For Two and Five phase motors, these 2 registers set the maximum step rate. If this value is set too high the motor may skip steps.

For DC Encoded motors, these 2 registers set the maximum torque.

WARNING: Changing these values from the factory defaults may result in a system that does not perform correctly.

See the document **Motor Controller User's Manual-Windows** for the proper factory defaults.

8.2.4.REG_SETUP_REVBACKLASH

REG_SETUP_REVBACKLASH_1(0x18) sets motor 1 reverse backlash value

REG_SETUP_REVBACKLASH_2(0x28) sets motor 2 reverse backlash value

These 2 registers set the number of position count to use for backlash when traveling in the reverse direction. This value will have no effect if the motor is setup to reverse seek thru home.

8.2.5.REG_SETUP_FWDBACKLASH

REG_SETUP_FWDBACKLASH_1(0x19) sets motor 1 forward backlash value

REG_SETUP_FWDBACKLASH_2(0x29) sets motor 2 forward backlash value

These 2 registers set the number of position count to use for backlash when traveling in the forward direction. This value is normally set to zero.

8.2.6.REG_SETUP_CONFIG

REG_SETUP_CONFIG_1(0x1B)	sets motor 1 sensor configuration
REG_SETUP_CONFIG_2(0x2B)	sets motor 2 sensor configuration

These 2 registers setup the sensor configuration and homing operation.

- Bit 0: 0 – near sensor is home
1 – far sensor is home
- Bit 1: 0 – reverse seek direct
1 – reverse seek thru home
- Bit 2: 0 – Axis 2 enabled (only for AXIS 2)
1 – Axis 2 disabled (only for AXIS 2)

The default is for the near sensor to be home and to reverse seek thru home.

8.2.7.REG_SETUP_WRITE

REG_SETUP_WRITE(0x0E)	commits system register changes
-----------------------	---------------------------------

Writing a 0x1 to this register commits any and all changes made to the setup registers to flash memory. This allows the changes to survive beyond a power down in the controller.

Example:

```
// set near sensor to home and turn off reverse seek thru home
v = 0;
ConnectionWrite(REG_SETUP_CONFIG_1, &v)
v = 0;
ConnectionWrite(REG_SETUP_CONFIG_2, &v)

// commits the changes to flash memory
// so it survives the next controller power cycle.
v = 0x00000001;
ConnectionWrite(REG_SETUP_WRITE, &v);
```

8.3. Motor Controller Operation Registers

These user registers are intended for normal operation.

8.3.1.REG_USER_TARGET

Writing to either of these two registers commands the motor to seek to the specified value.

REG_USER_TARGET_1(0x10) motor 1 seek command

REG_USER_TARGET_2(0x20) motor 2 seek command

Example: See the appropriate code sample.

8.3.2.REG_USER_INCREMENT

Writing to either of these two registers commands the motor to advance to retreat by the specified value. This command is intended for use to make minor adjusts to the position and do not use the backlash parameters or the seek thru home capabilities.

REG_USER_INCREMENT_1(0x11) motor 1 delta seek from current

REG_USER_INCREMENT_2(0x21) motor 2 delta seek from current

Example: See the appropriate code sample.

8.3.3.REG_USER_CURRENT

These two registers are used to query the current motor position.

REG_USER_CURRENT_1(0x12) motor 1 current position

REG_USER_CURRENT_2(0x22) motor 2 current position

Example: See the appropriate code sample.

Developers Interface Guide

8.3.4.REG_USER_LIMIT

Writes to these registers command the motors to do one of the following three operations:

REG_USER_LIMIT_1(0x13) motor 1 limit seek command

REG_USER_LIMIT_2(0x23) motor 2 limit seek command

The value written drives the motor to one of the limits:

- 0 – drive motor to home.
- 1 – drive motor to limit.
- 2 – abort current motor motion.

Example: See the appropriate code sample.

8.3.5.REG_USER_STATUS

Reads from these registers query the motor status.

The low 8 bits indicate the motion state. While the motor is in motion the lower 8 bits will return a non zero value.

REG_USER_STATUS_1(0x14)

REG_USER_STATUS_2(0x24)

Return values:

Lower 8 bits (Bit 0 thru 7).

Hex	Binary	Motor motion state
0x0	0000	0 – idle
0x1	0001	1 – driving to home
0x2	0010	2 – coming off home
0x3	0011	3 – driving to limit
0x4	0100	4 – seeking forward
0x5	0101	5 – deaccel in forward direction
0x6	0110	6 – forward backlash (i.e. rev after overshoot)
0x7	0111	7 – seeking reverse
0x8	1000	8 – deaccel in reverse direction
0x9	1001	9 – reverse backlash (i.e. forward after overshoot)
0xB	1011	11- forward deaccel during an abort
0xC	1100	12- reverse deaccel during an abort

Bit 8 – motor on home sensor

Bit 9 – motor on limit sensor

When testing for the end of motion, just the lower 8 bits should be compared to zero. If the entire word is compared zero the motion may be completed but we could be stopped on the limit sensor. The return value in this case would be 0x0000000200.

Example: See the appropriate code sample

8.3.6.REG_SETUP_LIMIT

Reads from these registers return the values at the number of position counts between HOME and LIMIT.

REG_SETUP_LIMIT_1(0x1C) returns motor 1 limit value

REG_SETUP_LIMIT_2(0x2C) returns motor 2 limit value

This value is measured during the power up motor calibration sequence.

9. BrightLight Register Definitions

9.1. BrightLight Identification Registers

9.1.1.REG_SYS_PRODUCTID

REG_SYS_PRODUCTID(0x01) read only

Return value:

Either

DEF_BRIGHTLIGHT 0x4000

DEF_MOTOR_CONTROLLER 0x4001

9.1.2.REG_SYS_VERSIONHW

REG_SYS_VERSIONHW(0x02) read only

Return value:

32 bit hardware version identifier.

9.1.3.REG_SYS_VERSIONDATE

REG_SYS_VERSIONDATE(0x03) read only

Return value:

32 bit build date identifier.

0x00,YY,MM,DD

Example:

0x00050402 would be 2005, April, 2

9.1.4.REG_SYS_VERSIONSW

REG_SYS_VERSIONSW(0x04) read only

Return value:

32 bit software version identifier.

9.2. BrightLight Operation Registers

9.2.1.REG_PWM_ABSOLUTE

REG_PWM_ABSOLUTE(0x40) read/write

This value is in 0.1 increments

0 -> 0%

1000 -> 100%

9.2.2.REG_PWM_INCREMENT

REG_PWM_INCREMENT(0x41)

write only

+/- change to current PWM value

This value is in 0.1 increments

10 -> 1%

-10 -> -1%